

Loan Library Catalogue

No-Tech Communication



E-Tran Frame

Designed for communication through eye pointing.

Clear Perspex frame with central hole to allow the communication partner to see where the communicator is looking.



PECS Book

Picture Exchange Communication System book with removable pages and sentence strip.

Add your own symbols/photos.

Low-Tech Communication



BIGmack

Single message recordable switch.

Record up to 2 minutes.

Large activation surface for playback.



Talking Tins

Single message recordable switch.

Record up to 20 seconds.

Available with adaptors and straps.



Big Step-by-Step

Multiple message recordable switch.

Record up three different sets of sequential messages.

Total of 4 minutes recording time.

Large activation surface for playback.



Little Step-by-Step

Multiple message recordable switch.

Record up three different sets of sequential messages.

Total of 4 minutes recording time.

Smaller activation surface for playback.



Adaptation Randomizer

Multiple message recordable switch.

Different modes of operation: single message, jumble (plays messages in random order) and eliminator (plays messages in random order but will not repeat message until all have been played)

Switches



Big Buddy Button

Large activation surface area (11.5cm).

Push to activate.

Click on activation.

Wired switch.



Jelly Bean

Small activation surface area (6.5cm).

Push to activate

Click on activation.

Wired switch.



Specs Switch

Very small activation surface area (3.5cm).

Push to activate.

Click on activation.

Wired switch.



Wobble Switch

Activated by fine or gross motor movements.

Wobble in any direction to activate.

Can be used by any part of the body including the head.

Wired switch.



Grasp Switch

Activate by holding in the palm of the hand and squeezing the rubber hand grip.

Robust and durable.



Candy Corn

Activated by proximity, needs no pressure to activate.

Medium sized surface is suitable for most users.

Give auditory and visual feedback when activated.

Environmental Controls



Switch Latch and Timer

Flexible appliance controller.

Multiple operation modes.

Works with most wired switches.



Jelly Beamer

Wireless transmitter and receiver.

Beamers pair with own receiver, so multiple can be used in the same room.

Effective range up to 10 meters.

Mounts



Flexzi 1 Mount

Heavy duty clamp for switches and small devices.

Single arm, adjustable length.

Clamp base will attach to a range of surfaces, including wheelchairs.




Flexzi 2 Mount

Heavy duty clamp for switches, small and heavier devices.

Double arm, adjustable length.

Clamp base will attach to a range of surfaces, including wheelchairs.



Maxess switch tray

Cushioned tray.

Will hold devices and switches in almost any position.

Different sizes.



Maxess Switch Mounts

Offers variety of positions for switches.

Right angle or sloping.

Velcro for easy mounting.